





 vacuum-based, hollowtubed instrument with a guarded, oscillating blade



 self-contained electrical machine; no grounding pad is needed

Microdebrider Advantages

- Removes both blood and tissue
- Able to be used by a single surgeon
- Small footprint in the operating room
- Absence of thermal injury
- Cost-effective







Microcebrider <u>Disadvantages</u>

- Unable to provide hemostasis
- Complications may occur very quickly
- The weight and vibration of the powered instrumentation may hinder tactile feedback
- Microdebrider adds unique complications

Microdebrider Rhinology Complications

A survey presented at AAO-HNS in 2004 showed...

- 1/3 of participants were unaware of orbital injury
- 1/2 of significant bony complications

...that may occur with powered instrumentation during sinonasal surgery.









There is a paucity of literature describing microdebrider complications in Laryngologic and Airway surgery....

Electrical Complications

- Intraoperative EKG waveform of 30y/o M resembled that of ventricular tachycardia; pt clinically stable...
- Intraoperative development of ventricular tachycardia in 25y/o F...



- Mortensen and Woo describe anterior glottic web following excision of a vocal fold polyp with a microdebrider
- Occasionally seen in several voice centers...%?
- See decreased mucosal wave

propagation/adynamic segment on stroboscopy



- U/L 40% PTX, pneumonediastinum & retroperitoneal air after papillomas removal with microdebrider and jet ventilation.
- Tri-cut microdebrider blade caused small integrity breaches in the trachea worsened with transtracheal jet ventilation.





- Nine Laryngology Centers
- >140 years patient care experience
- Discussion of major complications involving microdebrider

• Moderate hemorrhage, scarring of lamina propria and moderate webbing excluded



- 4 Centers had no major complications
- Inadvertent tissue removal
- Airway perforation
- Severe hemorrhage







- Subglottic anterior airway tear
- Tracheoesophageal fistula
- Tracheal tear, tx thoractomy/local muscle flap



 Three cases of severe hemorrhage requiring packing and angiography

One led to embolization led to CVA and death

 Intraop innominate artery violation resulted in death



- Pick the correct instrument
- Any concern for malfunction stop!
- Maintain clear and unobstructed view
- Judicious hemostasis
- Avoid high suction
- Be aware of this possibility

**Every major complication involved the Tricut Blade





